

# INDIAN INSTITUTE OF TECHNOLOGY

DATE 30/01/25

Expt: 6 Speed Viscometer (FANN-35S)

SHEET NO. 15

Aim: To determine the viscosity of a given fluid sample using 6 speed viscometer.

## Procedure

→ Ensure that the bob

Theory: FANN-35S 6-speed viscometer is used to determine the viscosity of a given fluid sample by measuring resistance to rotational motion at different speeds (RPM). Viscosity is a key property in fluid mechanics, particularly in drilling fluid analysis, where it impacts flow behaviour and performance.

## Working Principle

The rotor rotates inside the fluid sample, and the torque on the ~~bob~~ bob (attached to the bob shaft) is measured.

The dial readings at different speeds (600 RPM, 300 RPM, etc.) help determine viscosity characteristics using standard eqns.

→ Apparent viscosity (AV): 600 RPM reading / 2

→ Plastic Viscosity (PV): (600 RPM reading) - (300 RPM reading)

→ Yield Point (YP) = (300 RPM reading) - (PV)

## INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 16

### Procedure


- 1) Ensure the bob and sleeve is attached with the bob shaft and rotor.
- 2) Connect the instrument to a power source and turn it on.
- 3) Plug the FANN 35 motor to 110 volt voltage transformer
- 4) Fill the sample cup to 75% capacity (170ml)
- 5) If temp. control is required, use a thermocup and thermo-meter.
- 6) Put the sample cup in the base of the instrument and raise the base till the sample touch the given marker.
- 7) Adjust RPM control as required.
- 8) Wait for stabilizing the reading and note down viscosity in CP with corresponding temp. and RPM.
- 9) Turn off equipment and disconnect it from the power source.
- 10) Remove and clean the thermocup, bob, bob shaft, sleeve and sample cup with water and tissue paper.

P.R.E.

Observation

@  $T = 26^\circ C$

$\theta(\text{deg})$	RPM	Dial Reading	$\tau(\text{Pa})$	$M(\text{Pas})$
1021.8	600	98	50.078	0.05
510.9	300	71	36.281	0.07
340.6	200	60	30.66	0.09
170.3	100	45	22.995	0.135
10.218	6	10	9.709	0.95
5.109	3	17	8.687	1.7



# INDIAN INSTITUTE OF TECHNOLOGY

DATE \_\_\_\_\_

SHEET NO. 17

## Calculation and Lab Exercise

→ Derive formula for Yield Point:

$$Y.P \left( \frac{\text{Pa}}{100 \text{ ft}^2} \right) = (300 \cdot \text{RPM})_{DR} - \text{Plastic Viscosity}$$

→ from the plot,

$$\text{slope of the line} = m = \frac{\theta_{600} - \theta_{300}}{11.703} = PV \text{ (in CP)}$$

eq' of line:

$$\frac{y - \theta_{300}}{\theta_{600} - \theta_{300}} = \frac{x - 300}{600 - 300}$$

$$\therefore y - \theta_{300} = \frac{PV(x - 300)}{300}$$

$$\therefore Y.P = \theta_{300} = \frac{PV(0 - 300)}{300} = -PV$$

$$\therefore \boxed{Y.P = \theta_{300} - PV} \quad (\text{Proved})$$

$$\rightarrow \tau(\text{Pa}) = D.R \times 0.511 \quad \begin{matrix} \downarrow \\ \text{Spring} \end{matrix} \quad \begin{matrix} \dot{\gamma} \\ (\text{sec}^{-1}) \end{matrix} = \text{RPM} \times 1.703$$

$$\mu(\text{Pa.s}) = \frac{\tau}{\dot{\gamma}}$$

sample calculation for obs 1: RPM = 600, DR = 98

$$\therefore \boxed{\tau = 50.078 \text{ Pa}}$$

$$\dot{\gamma} = 600 \times 1.703 = \boxed{1021.8 \text{ sec}^{-1}} \quad M_2 \frac{50.078}{1021.8} = \boxed{0.05 \text{ Pa.s}}$$

$$\boxed{AP = \frac{98}{2} = 49 \text{ CP}}$$

$$\boxed{PV = 08 - 71 = 27 \text{ CP}} \quad \boxed{Y.P = 44 \frac{\text{Pa}}{100 \text{ ft}^2}}$$

P.R.E.

## INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 18

### Results

Apparent viscosity = 49 CP      Yield Point = 44 lb/100 ft<sup>2</sup>  
Plastic viscosity = 27 CP

### Precautions

### Discussion

- 1) Always ensure that the fluid temp. doesn't go beyond 180°F
- 2) Don't move "Gear Shift knob" if "viscometer switch" is off.
- 3) Don't try to rotate the Thermo-WP knob beyond marked position.
- 4) The rotor, bob and splash guard should be cleaned and dried after each operation.

### Conclusion

- i) The Fann-35S viscometer accurately determines fluid viscosity at varying rates.
- ii) Plastic viscosity indicates particle interaction and friction, while Yield Point reflects the fluid's ability to suspend solids.
- iii) This experiment highlights the importance of proper temp. control, equipment calibration, and fluid sample preparation to ensure accurate viscosity readings.
- iv) Future improvements could include automated data logging for efficiency and minimizing ~~for~~ manual reading error.

P.R.E.

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO

## Assignment

- 1) Temp.      2) To Protect bob shaft      3) Rotor

$$4) P.V. = (\text{DR})_{600} - (\text{DR})_{300} = 36 - 12 = 24 \text{ Pa}$$

$$\gamma P = (\text{DR})_{300} - PV = 12 - 24 = -12 \text{ lb/inch}^2$$