

# INDIAN INSTITUTE OF TECHNOLOGY

DATE 16/01/25

Exp 1: Introduction; Design and Preparation of WBM

SHEET NO. 1

Aim: Design and Preparation of WBM drilling fluids

Definition: Drilling fluid also called drilling mud, used to aid the drilling of oil and natural gas wells and on exploration drilling rigs.

Drilling fluid composition: oil-Based muds (<5% water),  
water-Based muds, Foams, Gas fluids, Reverse emulsions  
(5 to 50% water)

Drilling fluid additives:

⇒ Water Base:

- Bare fluid: clear water or brines
- Viscosifiers: clay (bentonite, hectorite), Polymers
- Filtration control agents
- Shale stabilizers
- weighting agents: low gravity solids (~2.7 sg), high gravity solids (~4.2 sg)

⇒ Oil Base:

- Bare fluid: isomeric olefins, diesel, esters etc.
- Emulsifiers
- water (or Brine)
- Viscosifiers: organophilic clays, Polymers
- Filtration control agents
- weighting agents: low gravity solids, High gravity solids

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 2

Equipment used: Hamilton Beach mixture

Standard operating procedure of Hamilton Beach Mixture:

- i) Formulate the required amount of water, salt, viscosifier, weighting material, Additives and LGS as req.
- ii) Clean/wipe the mechanical stinner properly before inserting into mixing container.
- iii) Measure all liquid components using measuring cylinder and solid components using properly balanced weighting balance for accuracy.
- iv) Add the ~~reqd~~ amount of water into mixing container.
- v) Add the ~~reqd~~ amount of salt into water and mix @ low speed for 90 sec.
- vi) Add the ~~reqd~~ amount of Viscosifier, Additives, LGS and weighting material in proper order with a time interval of 5 min after every addition of single component.
- vii) Continue the mixing until uniformity is obtained
- viii) Once uniformity is obtained stop mixing and remove the contents to separate container
- ix) Rinse and clean the mechanical stinner properly.
- x) Make sure RPM kept to zero before turning the mechanical stinner OFF and then power OFF main switch.

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 3

## Calculation and Results

Given data,

- Water based Mud  $\rightarrow P_f = 11 \text{ lb-bbl} (350 \text{ ml})$
- Fluid density  $\leftrightarrow 11 \text{ ppg} = 1.32 \text{ g/ml}$
- Salt = 7500 ppm = ~~75 g/L~~ = 7.5%
- Barite density = 4.2 SG

Formula used:

$$\begin{aligned} & \frac{P_f}{8.345} = V_{\text{Brine}} P_{\text{Brine}} + V_{\text{Add}} P_{\text{Add}} + V_{\text{Barite}} P_{\text{Barite}} \quad \left| \frac{V_{\text{Br}}}{V_w} = \frac{100 / P_{\text{Br}}}{(100 - \text{salt-wt.})} \right. \\ & V_f = V_{\text{Br}} + V_{\text{Barite}} + V_{\text{Add}}. \\ & M_{\text{Br}} = M_w + M_{\text{salt}} \end{aligned}$$

From additive details:

$$V_{\text{Add}} = \frac{1}{1} + \frac{2}{1.2} + \frac{15}{2.6} = 8.4358 \text{ ml/350 ml.}$$

$$\begin{aligned} & \text{from salt table, for NaCl} \quad P_{\text{Add}} = \frac{1+2+15}{(8.4358/350)} \\ & \text{for weight } \% = 7\%. \quad \text{Density} = 1.045 \text{ kg/L} \\ & \text{weight } \% = 8\%. \quad " = 1.056 \text{ kg/L} \end{aligned}$$

$$\therefore P_{\text{Brine}} = 1.052 \text{ kg/L for } 7.5\%.$$

$$\frac{V_{\text{Br}}}{V_w} = \frac{100 / 1.052}{(100 - 7.5)} = 1.0276$$

$$V_f = 350 = V_{\text{Brine}} + 8.4358 + V_{\text{Barite}}$$

$$\therefore V_{\text{Brine}} + V_{\text{Barite}} = 341.56 \quad (i)$$

$$V_{\text{Brine}}^f + V_{\text{Barite}}^f = 0.0758$$

$$\text{Now, } V_{\text{Brine}} \times 1.052 + \frac{8.4358 \times 15+2+1}{8.4358} + V_{\text{wt}} \times 4.2 = \frac{10}{8.345} 350 \times 1.3$$

$$\Rightarrow 1.052 V_{\text{Br}}^f + 4.2 V_{\text{wt}}^f = 1.147 \Rightarrow 1.052 (0.0758 - V_{\text{Br}}^f) + 4.2 V_{\text{wt}}^f$$

$$\therefore 462 = 1.052 V_{\text{Br}} + 4.2 V_{\text{Br}} + 18 \quad (ii)$$

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 4

$$\begin{aligned} V_{wt}^f &= \frac{0.1204}{3.148} = 0.038 \quad \text{On Solving (i), (ii) ---} \\ V_{Br}^f &= 0.0379 \quad \frac{V_{Br}}{V_w} = 1.0276 \Rightarrow V_w = 0.013 \\ V_{wt} &= \frac{26.899 \text{ ml}}{13.3 \text{ ml}} \quad V_{Br} = \frac{314.66 \text{ ml}}{318.66 \text{ ml}} \quad V_{add} = 8.4358 \text{ ml} \end{aligned}$$

Mass of Barite  $\rightarrow \cancel{13.3 \times 4.2 / 55.86} \quad 112.976 \text{ g.}$

$$\begin{aligned} \text{Mass of water} &\rightarrow \cancel{\rho_{\text{Brine}} \times 350 \times P_w - 0.013 \times 350 \times 1} \\ &\rightarrow \frac{314}{1.0276 \times 350} \times 350 = 305.56 \text{ g} = 319.45 \text{ g} \end{aligned}$$

$$\begin{aligned} \text{Salt to be added} &= M_{\text{Brine}} - M_{\text{water}} = 331.029 - 305.56 \\ &= \cancel{\rho_{\text{Brine}} (350 \times V_{Br}^f)} - 319.45 \\ &= 25.44 \text{ g} \end{aligned}$$

Results:

- 1) Mass of Barite: 112.976 g
- 2) Mass of NaCl: 25.44 g
- 3) Mass of Water: 305.56 g

## Precautions:

- 1) Continue to mix until uniformity is obtained in the solution
- 2) Rinse and clean the mechanical stirrer properly
- 3) Make sure to turn off the stirrer after use.

# INDIAN INSTITUTE OF TECHNOLOGY

DATE 16/01/25

Exp: 2 | Mud Balance

SHEET NO. 5


## Working Principle of Mud balance:

- A cup of known volume and man is attached to one end of the beam.
- Engineers fill the cup with drilling mud.
- The man of the filled cup is balanced on the other end of the beam by a fixed counterweight and a rider that can move freely along the graduated scale.



$$T_1 = F_1 \cdot a, T_2 = F_2 \cdot b$$

$$M_1 \cdot a = M_2 \cdot b$$



$F_1$  = input force to the lever

$F_2$  = output force

$a, b$  = is distance between the forces and the fulcrum.

Theory: Density is defined as weight per unit volume. It is expressed either in pounds per gallon (lb/gal) or pounds per cubic foot ( $\text{lb}/\text{ft}^3$ ), or in kg per cubic meter, or as specific gravity (SG). Pressure exerted by a static mud column depends on both the density and the depth. Mud density is used to control subsurface pressures and stabilize the wall bore.

## Procedure:

- 1 Remove lid from the cup, and completely fill the cup with the mud to be tested
- 2 Put the lid on the top of the cup and rotate until firmly seated, making sure mud is expelled through the hole in the lid.
- 3 Wipe the mud from the outside of cup by using clean tissue paper

DATE

4. Plan  
knife

5. Move  
by the  
be in

6. At the  
density

Precaution

- \* Use se
- \* Clean e

Vibration

checked w  
give a reo  
8.34 on th  
b the mud  
ater then  
new lid  
calibr

Calibration

Fusion

mud balance  
drilling mud,  
olling forma  
stance of Pr  
itegity and  
It! The




Fig : Mud Balance device

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 6

4. Place the balance arm on the base, with the knife-edge resting on the fulcrum.
5. Move the slide until the graduated arm is level, indicated by the level bubble present on the beam. (Bubble should be in centre of the markings and steady).
6. At the left-hand edge of the rider, note down the density value with temperature without disturbing the rider.

## Precautions:

- 1) Use separate dustbin to throw mud.
- 2) Clean everything after use.

Calibration Procedure: The mud balance calibration can be checked using fresh water. At (21C) fresh water should give a reading of 1.00 on the specific gravity scale psi/1000ft, 8.34 on the lb/gd scale, and 62.3 on the lb/in ft. scale.

→ if the Mud Balance doesn't give the correct reading for water then replacing the lid on the balance cup with a new lid can cause the Mud Balance to be out of calibration.

Calibration frequency: once a month.

## Conclusion:

The mud balance experiment is essential for determining the density of drilling mud, which plays a crucial role in wellbore stability and controlling formation pressures. This experiment reinforces the importance of proper mud weight management for overall well integrity and operational success.

Result: The density of Mud was WBM was found to be P.R.E. 11 ppg from the given exp.

## Assignment 01 (A)

→  $V_w$  = vol. of water req.       $V_{Br}$  = vol. of Brine

$w_s$  = mass of salt in weight %.

$\rho_{Br}$  = density of brine

Brine's total weight ( $W_{Br}$ ) per unit vol.  $\Rightarrow W_{Br} = \rho_{Br} \times V_{Br}$

weight of salt :  $w_s = \text{salt-wt\%} \times W_{Br}$

$$\therefore w_s = \text{salt-wt\%} \times (\rho_{Br} \times V_{Br})$$

weight of water:  $W_w = (100 - \text{salt-wt\%}) \times W_{Br}$

$$W_w = (100 - \text{salt-wt\%}) \times (\rho_{Br} \times V_{Br})$$

~~Also~~,  $W_w = V_w \times 8.345$

$$\Rightarrow V_w \times 8.345 = (100 - \text{salt-wt\%}) \times (\rho_{Br} \times V_{Br})$$

Now,  $\frac{V_{Br}}{V_w} = \frac{1}{\rho_{Br}} \cdot \frac{100}{(100 - \text{salt-wt\%})} \times 8.345$

$$\therefore \boxed{\frac{V_{Br}}{V_w} = \frac{\frac{100}{\rho_{Br}/8.345}}{100 - \text{salt-wt\%}}}$$

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 8

2)

## Dispersed System

Definition: A system where fine solids are intentionally dispersed into the fluid using chemical dispersants.

Key Additives: Dispersants (e.g. lignosulfonates, tannins) to break down and evenly distribute particles.

Fluid characteristics: Higher viscosity due to dispersed fine solids

Solid control: Solids are dispersed to prevent sedimentation

Applications: Used for wells with high temp. and high pressure (HTHP) or in reactive clay formations.

## Non-dispersed System

System where fine solids are not dispersed but allowed to settle or be removed.

Minimal chemical additives, typically uses viscosifiers (bentonite)

lower viscosity, as fine solids are not dispersed

Solids are controlled via mechanical removal or settling.

Used in simple, low-cost wells with minimal reactive clays.

3) Arrangement of additives in typical order of addition for WBM formation:

Salt → ~~viscosifier~~ weighting Material → ~~viscosifier~~ weighting Material

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO.

9

## Assignment 01 (B)

- 1)  $N/m^3$  not a unit of density, unit of specific density.
- 2) Density of substance changes with temperature, so it is necessary to mention the temp at which specific gravity is calculated.
- 3) True; because it can lead to inconsistencies in the measurement.
- 4) Hydrostatic Pressure (Psi) = Mud weight (ppg)  $\times$  TVD (ft)  $\times$  0.052  
=  $11 \times 9000 \times 0.052 = 5148$  Psi

# INDIAN INSTITUTE OF TECHNOLOGY

DATE 23/01/25

EXP 3: Digital pH Meter

SHEET NO. 10

Aim: To determine the pH of a given fluid using pH meter.

## Theory:

pH of the solution is the logarithm of the reciprocal of the ( $H^+$ ) concentration in gm/lit.

The acidity and the alkalinity of the drilling fluid can be measured by the conc. of the ( $H^+$ ) ion in the fluid.

$$pH = -\log[H^+] = \log \frac{1}{[H^+]}$$

on the scale 7 is neutral, less than 7 is acidic and greater than 7 is alkaline.

## Procedure:

- 1) Before measurement, rinse pH electrode and ATC probe with distilled water to remove any impurities stuck onto the bodies of probes.
- 2) Wipe pH electrode and ATC probe with clean dry tissue paper
- 3) Power on the meter using ON/OFF key. Press 'MODE' key to select your desired mode (pH) of operation (pH, mV, Ion, or temperature)
- 4) Dip and stir both probes gently into given sample, swirl gently and wait for the reading to stabilize. Note the reading. Freeze the display if desired.
- 5) Note the value of pH with temp.
- 6) Rinse probe with distilled water before taking next reading or storage

DATE




fig: pH Measurement

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 11

## Working Principle of pH Meter:

- 1) The  $H^+$  ions Present
- 2) The working principle of pH meter is based on the relation between the electric voltage and hydrogen ion concentration.
- 3) The pH meter consists of glass (also called as measurement electrode) and reference electrode.
- 4) Both electrode are connected to a voltmeter
- 5) The measurement electrode consists of a glass membrane, which is sensitive to hydrogen ion concentration.
- 6) The reference electrode is standard and has constant potential.
- 7) The acidic solution is rich in  $H^+$  ion concentration
- 8) When pH probe is dipped in an acidic solution and switch on the voltmeter
- 9) The  $H^+$  ions move close to the glass membrane of the sensitive glass bulb (external side of the bulb)
- 10) A similar reaction occurs inside the bulb, which is filled with a buffer solution of neutral pH (which has constant no. of hydrogen ions).
- 11) Hence, this causes the difference in the concentration of hydrogen ions across the membrane
- 12) When hydrogen ion conc. inside bulb is less than outside solution (test solution), then given solution is acidic and  $pH < 7$ .

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 12

- 13) When hydrogen ion conc. inside glass bulb is similar to the outside solution, then given solution is natural pH hence the pH is 7.
- 14) When hydrogen ion concentration inside the glass bulb is more than the outside solution (test sol.), the given solution is alkaline and  $\text{pH} > 7$ .

## Calculation and Results

- 1) pH = 7.16 for Distilled water @  $24.5^{\circ}\text{C}$  ✓
- 2) pH = 7 for buffer solution (Given)  
observed pH = 7.04 @  $24.7^{\circ}\text{C}$  ✓
- 3) for pH = 10 buffer solution  
observed pH = 9.8 @  $24.8^{\circ}\text{C}$  ✓
- 4) for WBM, observed pH = 7.30 @  $24.7^{\circ}\text{C}$  ✓

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 13

## Precautions:

- 1) Don't use the electrode as a stirring glass rod when adjusting pH.
- 2) Do the calibrations before testing the each sample using standard buffer solutions.
- 3) Avoid temp. fluctuations and never keep the pH meter exposed to direct sunlight.
- 4) Always do the calibration by attaching ATC to the electrode for best accuracy.
- 5) Glass electrode shouldn't be left out of the storage solution for prolonged intervals as the glass membrane get dehydrated resulting in slower response or response failure.

## Discussion

The pH of the mud sample was found to be 7.30 at 24.7°C, indicating that it is slightly alkaline. This suggests that the basic components in the sample, making it more neutral but leaning towards alkalinity. Proper pH balance in drilling mud is ~~an~~ important for maintaining the drilling efficiency and preventing corrosion or sedimentation of the equipment.

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 14

## Assignment

- 1) Full form of ATC in the pH measurement equipment:  
Automatic Temp. Compensation ✓
- 2) pH meter shouldn't be exposed to direct sunlight because:  
Effect on accuracy of the equipment. ✓
- 3) During calibration of pH meter buffer sol. of pH 7 should be taken first.
- 4)  $\text{pH} = -\log[\text{H}^+]$  | given:  $[\text{H}^+] = 10.45 \times 10^{-3} \text{ M}$

$$\begin{aligned}\therefore \text{pH} &= -\log(10.45 \times 10^{-3}) = -\log(10.45) + 3 \\ &= -1.02 + 3 = 1.98\end{aligned}$$

# INDIAN INSTITUTE OF TECHNOLOGY

DATE 30/01/25

Expt: 6 Speed Viscometer (FANN-35S)

SHEET NO. 15

Aim: To determine the viscosity of a given fluid sample using 6 speed viscometer.

## Procedure

→ Ensure that the ~~bob~~

Theory: FANN-35S 6-speed viscometer is used to determine the viscosity of a given fluid sample by measuring resistance to rotational motion at different speeds (RPM). Viscosity is a key property in fluid mechanics, particularly in drilling fluid analysis, where it impacts flow behaviour and performance.

## Working Principle

The rotor rotates inside the fluid sample, and the torque on the ~~rotor~~ bob (attached to the bob shaft) is measured.

The dial readings at different speeds (600 RPM, 300 RPM, etc.) help determine viscosity characteristics using standard eqns.  
→ Apparent viscosity ( $\eta_v$ ):  $600 \text{ RPM reading} / 2$   
→ Plastic viscosity ( $\eta_p$ ):  $(600 \text{ RPM reading}) - (300 \text{ RPM reading})$   
→ Yield Point ( $\eta_y$ ) =  $(300 \text{ RPM reading}) - (\eta_p)$






fig: Schematic of FANN-35S Viscometer

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 16

## Procedure


- 1) Ensure the bob and sleeve is attached with the bob shaft and rotor.
- 2) Connect the instrument to a power source and turn it on.
- 3) Plug the FANN 35 motor to 110 volt voltage transformer
- 4) Fill the sample cup to 75% capacity (170ml)
- 5) If temp. control is required, use a thermocup and thermo-meter
- 6) Put the sample cup in the base of the instrument and raise the base till the sample touch the given marker
- 7) Adjust RPM control as required.
- 8) Wait for stabilizing the reading and note down viscosity in CP with corresponding temp. and RPM.
- 9) Turn off equipment and disconnect it from the power source.
- 10) Remove and clean the thermocup, bob, bob shaft, sleeve and sample cup with water and tissue paper.



Observation

@  $T = 26^\circ C$

$\delta (\text{Pa})$	RPM	Dial Reading	$\tau (\text{Pa})$	$\mu (\text{Pa.s})$
1021.8	600	98	50.078	0.05
510.9	300	71	36.281	0.07
340.6	200	60	30.66	0.09
170.3	100	45	22.995	0.135
10.218	6	10	9.709	0.95
5.109	3	17	8.687	1.7



$$\mu (\text{Pa.s})$$

$$\text{sample C}$$

$$\delta = 6$$

$$\frac{\theta}{2} = 4$$

# INDIAN INSTITUTE OF TECHNOLOGY

DATE \_\_\_\_\_

SHEET NO. 17

## Calculation and Lab Exercise

→ Derive formula for Yield Point:

$$Y.P \left( \frac{lb}{100 \text{ ft}^2} \right) = (300 \cdot \text{RPM})_{\text{OR}} - \text{Plastic Viscosity}$$

→ from the plot,

$$\cancel{\text{slope of the line} = m} \quad \frac{\theta_{600} - \theta_{300}}{7/300} = P.V \text{ (in CP)}$$

eqn of line:

$$\frac{y - \theta_{300}}{\theta_{600} - \theta_{300}} = \frac{x - 300}{600 - 300}$$

$$\therefore y - \theta_{300} = \frac{P.V(x - 300)}{300}$$

$$\therefore Y.P - \theta_{300} = \frac{P.V(0 - 300)}{300} = - P.V$$

$$\therefore Y.P = \theta_{300} - P.V \quad (\text{Proved})$$

$$\rightarrow \tau(Pa) = D.R \times 0.511$$

↳ Spring constant

$$\left( \frac{\text{N}}{\text{m}} \right) = \text{RPM} \times 1.703$$

$$M(\text{Pa.s}) = \frac{\tau}{\gamma}$$

Sample calculation for obs 1: RPM = 600, DR = 98

$$\therefore \tau = 50.078 \text{ Pa}$$

$$\gamma = 600 \times 1.703 = 1021.8 \text{ sec}^{-1} \quad M_2 = \frac{50.078}{1021.8} = 0.05 \text{ Pa.s}$$

$$AV = \frac{98}{2} = 49 \text{ CP}$$

$$P.V = 98 - 71 = 27 \text{ CP} \quad Y.P = 44 \frac{\text{lb}}{100 \text{ ft}^2}$$

### Results

Apparent viscosity = 49 CP      Yield Point = 44 lb/100 ft<sup>2</sup>  
 Plastic viscosity = 27 CP

### Precautions

### Discussion

- 1) Always ensure that the fluid temp. doesn't go beyond 180°F
- 2) Don't move "Gear Shift knob" if viscometer switch is off.
- 3) Don't try to rotate the Thermo-ATP knob beyond marked position.
- 4) The rotor, bob and splash guard should be cleansed and dried after each operation.

### Conclusion

- i) The Fann-35S viscometer accurately determines fluid viscosity at varying rates.
- ii) Plastic viscosity indicates particle interaction and friction, while Yield Point reflects the fluid's ability to suspend solids.
- iii) This experiment highlights the importance of proper temp. control, equipment calibration, and fluid sample preparation to ensure accurate viscosity readings.
- iv) Future improvements could include automated data logging for efficiency and minimizing ~~the~~ manual reading error.

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 19

## Assignment

- 1) Temp.      2) To Protect bob shaft      3) Rotor

$$4) P.V. = (DR)_{600} - (DR)_{300} = 36 - 12 = 24 \text{ Pa}$$

$$\gamma P = (DR)_{300} - PV = 12 - 24 = -12 \text{ lb/100 ft}^2$$

# INDIAN INSTITUTE OF TECHNOLOGY

DATE 05/02/25

EXP: 5 | Marsh funnel Viscometer

SHEET NO. 20

Aim: To measure the viscosity of drilling fluid using Marsh funnel Viscometer.

Theory: Marsh funnel Viscometer is a sample device used to measure the viscosity of drilling fluids. It operates on the principle that the time taken for a known volume of fluid to flow through a standardize funnel correlates with its viscosity.

The key theoretical aspects include:

1) Fluid flow: rate at which fluid flows through the funnel is inversely proportional to its viscosity.

2) Effective Viscosity: The viscosity measured is an effective viscosity, which is calculated using given formula:

$$\boxed{\mu = \rho(t-25)}$$

$\mu$  = effective viscosity in centipoise  
 $\rho$  = density in g/cm<sup>3</sup>  
 $t$  = quest funnel time in seconds

3) Temp. sensitivity: The viscosity is a function of temp, which is why temp. measurement is part of Procedure.

## Procedure:

- 1) Hold the funnel in an upright position just above measuring cup.
- 2) Block the orifice outlet with your index finger
- 3) Pour drilling mud through the screen until fluid level reaches.




Fig: Schematic of Marsh funnel.

$$M = 1.0$$

### Precautions

- Clean and use.
- Don't be

IND

DATE

- iv) Wait until working
- v) Remove and store
- vi) Measure the mass
- vii) Measure
- viii) Recom

### Calculations

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 21

- iv) Wait for 15 sec. and ensure your stopwatch is working properly.
- v) Remove your index finger from the orifice outlet and start the stopwatch simultaneously.
- vi) Measure the time it takes for drilling mud to fill the marked portion (946 ml or 1 quart) of measuring cup.
- vii) Measure temp. of fluid in  $^{\circ}\text{C}$  or  $^{\circ}\text{F}$ .
- viii) Record the viscosity in seconds and temp.

## Calculations and Result

$$\mu = P(t - 25)$$

$$\mu = 1.44 \times (44.5 - 25)$$

$$= 28.08 \frac{\text{g.s}}{\text{c.c.}} = 2808 \text{ cP}$$

## Observation

$$P_{212} \text{ PPG} = 1.44 \text{ g/c.c}$$

$$t = 44.5 \text{ sec}$$

$$T = 25^{\circ}\text{C}$$

## Precautions

- v) Clean and dry the funnel and cup thoroughly after each use.
- vi) Don't bend or flatten the orifice in the bottom.

Assignment

$$\Rightarrow 1500 \text{ cm}^3$$

$$\Rightarrow \text{Mesh size} = 12 \text{ inch}$$

~~3> Soln not~~  
3> All of the above

$$\Rightarrow M = \rho(t - 25) \quad \left| \begin{array}{l} \rho = \text{density} \\ t_2 = \text{quart funnel time} \end{array} \right.$$

$$\Rightarrow M_1 = 1.4(135 - 25) = 154 \frac{\text{g}}{\text{cm}^3 \cdot \text{sec.}}$$

$$\Rightarrow M_2 = 2.8 \times (110 - 25) = 238 \frac{\text{g}}{\text{cm}^3 \cdot \text{sec.}}$$

$$\Rightarrow M_3 = 2.5 \times (28 - 25) = 7.5 \frac{\text{g}}{\text{cm}^3 \cdot \text{sec.}}$$

DATE 13/02/25

Ex6: LPLT filter Press

SHEET NO. 23

Aim: To determine the fluid viscosity loss using ~~LPLT~~ <sup>LPLT</sup> API filtration apparatus.

Theory: The LPLT filter press is used to measure the filtration characteristics of drilling fluids under low temperature and low pressure conditions. The test simulates the formation of filter cake on the wellbore wall and the invasion of mud filtrate into the formation. In this process while applying a pressure of 100 psi to a mud sample for 30 mins, during which the liquid phase (filtrate) passes through a standardized filter paper. The volume of the filtrate collected and the thickness of the filter cake formed are key indicators of the mud's ability to control the fluid loss and maintain wellbore stability.

Procedure:

- 1) Prepare the sample according to API specifications. Stir the test sample atleast 10 mins with a high speed mixer.
- 2) Assemble the filtration cell first to fit the neoprene gasket into the base cap.
- 3) Fill the cell with approximately 400 ml of mud sample, leaving 1-1.5 cm of void space.
- 4) Place the assembly on the frame.
- 5) Install a CO<sub>2</sub> cartridge in the barrel and tighten it.
- 6) Tighten T-screw to fix assembly in the frame.
- 7) Position a dry graduated cylinder under the filter tube.

DATE

IND.

- 8) Ensure the
- 9) Run the
- 10) Carefully
- Take this

### Observation

	time (min)
	7.5
	10
	20
	30

### Result

... Total flow

The total

We have C  
eg of time:

150, the mud

### Cautions


Before beginning  
and dry, Par

Be Careful

Maximum  
beyond of  
Should not

Fig: Top view of Well Bore

(Mud Filtrate Invasion)



- 8) Ensure the bleeder valve is closed.
- 9) Run the exp. for 30 mins, measuring the filtrate volume.
- 10) Carefully remove the filter paper and measure the filter cake thickness in mm.

### Observation

time(min)	Vol. (ml)
7.5	1.7
10	1.9
20	2.7
30	3.5

### Result

~~Total fluid loss = 3.5 ml~~

The total fluid loss has been shown in the above table.  
We have calculated the last two observation by extrapolation  
~~of time~~

Also, the mud cake thickness = 1mm.

### Precautions

- 1) Before beginning a test, make sure each part of the cell is clean and dry, particularly the screen.
- 2) Be careful while adjusting the pressure
- 3) Maximum allowable pressure is 100 psi. Don't operate beyond this pressure.
- 4) Should not expose the carbon cartridges to heat.

# INDIAN INSTITUTE OF TECHNOLOGY

SHEET NO. 25

DATE

## Assignment

1) Residue: solid portion after filtration ✓

2) Filtration area = 7.1 inch<sup>2</sup> ✓

3) Correct sequence:

Neoprene gasket → Mesh Screen → Filter paper → Neoprene gasket

4) Law behind filtration: Darcy's law. ✓

5) Amount of filtrate obtained after 30 mins: 11 cm<sup>3</sup>

Thickness of filter cake: 2.5 mm

$$\text{Darcy's law: } \frac{dV_f}{dt} = \frac{KA\Delta P}{\mu h_{mc}}$$

$V_f$  = vol. of filtrate in cm<sup>3</sup>

t = time in second

K = average permeability of  
mud cake in Darcy

$$A = 45 \text{ cm}^2 / \Delta P = 100 \text{ psig} = 6.8 \text{ atm}$$

$$\text{Amount of filtrate obtained in 30 min} = 11 \text{ cm}^3$$

$$\text{Thickness of filter cake} = 2.5 \text{ mm.} \quad ?$$

$$K = \frac{11 \times 1.5 \times 2.5}{45 \times 30 \times 60 \times 6.8} = 7.47 \times 10^{-6} \text{ Darcy.} \quad | \begin{array}{l} \text{Darcy} \\ = 9.8 \times 10^{-13} \text{ m} \end{array}$$

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

Ex7: Total Hardness Determination  
and Chlorides Using Titration

SHEET NO. 26

Aim: Determination of Total Hardness (Ca) and Chlorides in drilling fluid using Titration Method

Theory: Hardness of water or mud filtrate is due to the primarily to the presence of Calcium and magnesium ions. The harder the water, the more difficult for many chemicals to function, particularly bentonite clay. Hardness decreases the efficiency of most polymers. Excessive hardness will cause flocculation of clays in the mud, effecting filtration rates.

Procedure:

→ Total Hardness Determination:

i) Add 25 ml deionized water to titration vessel, then add 1ml (approx. 10 drops) of Versenate Hardness Buffer solution.

ii) Add 3 drops of versenate Hardness Indicator solution. If calcium and/or magnesium are present in the deionized water, a wine red color will develop. Otherwise, will remain blue.

iii) While stirring, titrate with EDTA (Standard Versenate) until the solution color changes from wine-red to blue.

iv) Note ml of EDTA and calculate Hardness, mg/L

v) Add 1 ml filtrate to deionized water. If  $\text{Ca}^{2+}/\text{Mg}^{2+}$  ion are present wine red color will develop. While stirring, titrate with EDTA until color indicator changes to blue.

Observation  
~~Total Hardness of given filtrate is~~

~~Amount of EDTA sol. used for Total Hardness determination: 1.2 mL~~

$$\therefore \text{Total Hardness (calcium)} \text{ of given filtrate} = 1.2 \times 4000 \\ = 4800 \text{ mg/L}$$

DATE

NOW, TO FIND THE CHLORIDES PRESENT IN GIVEN DRILLING FLUID -

APPARATUS AND CHEMICAL / REAGENTS REQUIRED

- i) Burette      ii) Pipette      iii) Wash Bottle      iv) Spatula
- v) Conical flask      vi) Beaker      vii) Silver Nitrate Sol.
- viii) Potassium Chromate      ix) Mud filtrate sample

SIGNIFICANCE OF THE EXP:

- Chlorides (mg/L) is determined by the salinity of water-based mud.
- Chlorides act as a contaminant in freshwater mud system.
- Salt contamination can cause an increase in viscosity, gel strengths, and fluid loss.
- If you see changes in the chloride content, it can be indications of drilling in a saline formation, or to take the influx of water in reservoirs.

PROCEDURE

- 1) Add 1 ml of the filtrate to the titration vessel
- 2) Add 2-3 drops of phenolphthalein to the test sol.
- 3) If the pink color disappears add 0.2N sulphuric acid until pink color disappears
- 4) If no pink color is seen skip the  $H_2SO_4$  addition.
- 5) Add 5-10 drops of Potassium Chromate to titration vessel
- 6) while stirring, titrate with 0.282/0.1  $AgNO_3$  sol. till the first color disappears

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 28

## Observation

Amount of  $\text{AgNO}_3$  sol. used for titrating filtrate: 10 mL

## Calculation

if using 0.282 N  $\text{AgNO}_3$ :

$$\text{Chlorides (mg/L)} = 10 \times 10,000 = 10 \times 10^5 \text{ mg/L}$$

if using 0.1  $\text{AgNO}_3$ : chlorides (mg/L) =  $\frac{10 \times 10,000}{2.82} = 3546 \text{ mg/L}$

## Results

- Amount of Chlorides present in given filtrate sol = 3546 mg/L  
→ Total hardness of a given filtrate sol. = 4800 mg/L  
→ Amount of Calcium present in given filtrate = 4800 mg/L

## Precautions

- i) Make sure that there should not be any leakage from burette during titration.
- ii) Remove the air bubbles from the nozzle of the burette.
- iii) Use your index finger while pipetting the sol.
- iv) Keep your eye in level with the liquid surface while taking the burette reading or while reading the pipette or measuring flask etc.

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 29

Asgn: Total Hardness

- 1) The hardness of mud filtrate is due to the presence of Chlorides of Calcium and Magnesium.
- 2) Temporary hardness in water is due to the presence of Carbonates and Chlorides of calcium.
- 3) Permanent hardness of water is removed by all methods except Boiling Process.
- 4) Two Type of hardness: Temporary Hardness, Permanent Hardness.

Temporary hardness is caused by: Presence of dissolved bicarbonates and carbonates of calcium and magnesium.

Permanent hardness is caused by: Presence of sulfates, chlorides, or nitrates of  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$ .

- 5) Indicator used in determination of total hardness: Versenate Hardness Indicator.

DATE

### Asgn Chloride

- 1) Test is performed on Water Based Mud system
- 2) The chloride concentration is determined by titration with AgNO<sub>3</sub> solution.
- 3) End point of titration is detected by Potassium Chromate(K<sub>2</sub>CrO<sub>4</sub>) indicator.

### Asgn: Numerical

$$\text{Vol. of AgNO}_3 = 2.1 \text{ ml}$$

$$\text{filtrate vol.} = 1 \text{ ml}$$

$$\text{Normality of AgNO}_3 = 0.282 \text{ N}$$

$$\text{for chlorides (mg/L)}: 2.1 \times 10,000 = 21000 \text{ mg/L}$$

$$\text{Total Hardness} = \frac{\text{ml of EDTA} \times 4000}{\text{ml of filtrate}} = \frac{0.9 \times 4000}{1} = 3600 \text{ mg/L}$$

$$[\text{Ca}^{2+}] = \frac{40}{140} y = 3600 \Rightarrow y = \frac{3600 \times 110}{40} = 9900 \text{ ppm}$$

$$[\text{Cl}^-] = \frac{35}{58} x + \frac{70}{110} y = 21000$$

$$\Rightarrow \frac{35}{58} x + \frac{70}{110} \times 9900 = 21000 \Rightarrow x = 24360 \text{ ppm}$$

$$\therefore \text{NaCl amount} = 24360 \text{ ppm}$$

$$\text{CaCl}_2 \sim = 9900 \text{ ppm}$$

$$\left. \begin{array}{l} \text{Total Cl}^- = 21000 \text{ ppm} \\ \text{Ca}^{2+} = 3600 \text{ ppm} \end{array} \right\}$$

DATE 06/03/25

Exp : 8 | Mud filtrate Alkalinity Test

SHEET NO. 31

Aim: To find the alkalinity in the mud filtrate by titration method.

Theory: Alkalinity in drilling mud filtrate is primarily due to the presence of carbonate ( $\text{CO}_3^{2-}$ ), and bicarbonate ( $\text{HCO}_3^-$ ) ions. It plays a crucial role in maintaining pH, preventing corrosion, and ensuring stability of drilling fluids. The alkalinity measurement is typically conducted through acid-base titration, where the filtrate is titrated using phenolphthalein (P) and methyl orange (M) indicators to determine different forms of alkalinity.

Titration Process and Interpretation :


Phenolphthalein Alkalinity (P<sub>F</sub>): Measured when pH is above 8.5.

It indicates the presence of  $\text{CO}_3^{2-}$ .

Methyl Orange Alkalinity (M<sub>F</sub>): Measured when pH is between 4.3 and 8.5. It represents the presence of  $\text{CO}_3^{2-}$  and  $\text{HCO}_3^-$ .


Chemical Rxn involved :

1) Reaction of sodium carbonate with sulfuric Acid at higher pH (>8.5):



The change in color of phenolphthalein from pink to colorless indicates the conversion of all carbonate into bicarbonates.

2) Reaction of Sodium Bicarbonate with Sulfuric Acid (4.3 < pH < 8.5)



This conversion reaction is confirmed with the change in color of the filtrate in methylene blue titration test.

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 32

Apparatus : Pipette, burette, wash bottle, spatula, beaker, conical flask

Chemical/Reagents: mud, filtrate, distilled water, Phenolphthalein indicator, methyl orange indicator, sulphuric acid

## Procedure:

- Filtrate Alkalinity ( $P_F$ ):
- 1) Transfer 1 ml of LPLT filtration filtrate into the titration vessel using appropriate pipette.
  - 2) Add 3 drops of phenolphthalein indicator sol.
  - 3) If sol. turns pink, titrate (drop by drop) with 0.02 N sulphuric acid until pink color just disappears.
  - 4) Report the phenolphthalein alkalinity of filtrate  $P_F$ , as the number of  $\text{cm}^3$  of 0.02 N acid req. per  $\text{cm}^3$  of filtrate.
  - 5) Filtrate (if any)  $P_F$  alkalinity can indicate the presence of carbonate ions in most mud systems.
- Filtrate Alkalinity ( $M_F$ ):
- 1) To the same sample that has been titrated to the  $P_F$  endpoint, add 3 drops methyl orange indicator solution.
  - 2) Titrate with 0.02 N sulphuric acid (drop by drop) while stirring until the indicator changes from pale yellow to pinkish orange.
  - 3) Report the methyl orange Alkalinity,  $M_F$  as the total  $\text{cm}^3$  of 0.02 N acid per  $\text{cm}^3$  of filtrate req. to reach the methyl orange endpoint. The  $M_F$  is recorded as the total vol. of acid req. per  $\text{cm}^3$  of filtrate. (The total vol. indicates the  $\text{cm}^3$  of acid in the  $P_F$  titration plus the  $\text{cm}^3$  of acid  $M_F$  titration).

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 33

## Observation

for filtrate 1  
Amount of 0.02N sulphuric acid req. for filtrate

Alkalinity ( $P_F$ ): 0.8 ml. =  $T_1$

Amount of 0.02N sulphuric acid req. for filtrate Alkalinity ( $M_F$ ):

$$T_1 + T_2 = (0.8 + 0.7) \text{ ml} = 1.5 \text{ ml.}$$

## for filtrate 2

$$P_F = 0 = T_1$$

$$T_2 = 1.1$$

$$M_F = 1.1$$



## Calculation:

For, filtrate 1:  $P_F = 0.8$   $M_F = 1.5 \approx 2P_F$

$\therefore$  Carbonate ( $CO_3^{2-}$ ) =  $1200 \times P_F = 1200 \times 0.8 = 960 \text{ ppm}$   
only Carbonate is present.



For filtrate 2:  $P_F = 0$ ,  $M_F = 1.1$

$$\therefore [CO_3^{2-}] = 0 \quad [HCO_3^-] = 1220 \times M_F \\ = 1220 \times 1.1 = 1342 \text{ ppm}$$

## Precautions:

- Make sure there shouldn't be any leakage from burette.
- Remove air bubbles from nozzle of bottle burette.
- Use your index finger while pipetting the sol.

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 34

## Assignment

Ques: C

$$T_1 = 0.8 \text{ mL}$$

$$T_2 = 0.7 \text{ mL}$$

$$\therefore P_F = 0.8 \quad M_F = 0.8 + 0.7 = 1.5$$

$$\cancel{2} \cdot P_F = 2 \quad M_F \approx 2P_F$$

$\therefore [HCO_3^-] = 0$  only carbonate ( $CO_3^{2-}$ ) will be there.

$$3) [CO_3^{2-}] = 1200 \times P_F = 960 \text{ ppm}$$

$$4) [HCO_3^-] = 0 \text{ ppm.}$$

5) So, there is no bicarbonate present in the given filtrate solution. There is carbonate only, which can increase alkalinity, which may help control corrosion in drilling equipment but can also lead to scale formation in pipes and equipment.

K0  
10

# INDIAN INSTITUTE OF TECHNOLOGY

DATE 20/03/25

Exp 9: Cation Exchange Capacity

SHEET NO. 35

Aim: To determine the cation exchange capacity of a given drilling fluid sample using Methylene Blue Test.

## Theory:

- Methylene Blue capacity gives an estimate of the total cation exchange capacity (CEC) of the solids in the drilling fluid.
- Methylene Blue is a cationic dye, which strongly attracts to negatively charged sites on clays.
- The test, therefore, is a measure of the reactive clay concentration of a water-based drilling fluid.
- Increasing MBTs usually an indication of an increase in low gravity solids concentration.

## Procedure:

- Pipette 10ml of deionised water into the Erlenmeyer flask
- Use syringe to accurately measure 2ml of freshly agitated drilling fluid. Transfer this portion into flask.
- Boil this mixture on hotplate for 10 min. Don't boil to dryness.
- Add deionized water to dilute this mixture to approx. 50ml.
- Dip a stirring rod into the flask while the solids are suspended, and then remove it. Place a drop of the liquid on filter paper. The initial endpoint is reached when dye appears as a blue ring surrounding the dyed solid.
- When blue ring spreading from the spot appears, shake the flask an additional 2 mins. Place another drop on the filter paper. If the blue ring appears again, the final endpoint has been reached.

If blue ring doesn't appear, continue the titration until a drop taken after two minutes shows the blue tint. P.R.E.

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 36

## Observation

- Amount of Drilling fluid = 2 ml.  
→ Amount of Methylene Blue indicator = 2.5 ml.

## Calculation:

$$\rightarrow \text{CEC} = \frac{\text{Amount of Methylene Blue (ml)}}{\text{Amount of drilling fluid (ml)}} \\ = 2.5 / 2 = 1.25$$

$$\rightarrow \text{Bentonite, lb/bbl drilling fluid} = 5 * \text{Methylene Blue capacity} \\ = 5 * 2.5 = 12.5 \text{ ml.}$$

## Precautions

- 1) Make sure there is no leakage from burette during titration
- 2) Remove the air bubbles from nozzle of burette
- 3) Use your index finger while pipetting the solution.

## Assignment | Gen > C

- 1) Negatively Charged
- 2) Hydrogen peroxide used to remove organic materials, such as polymers, and thinning agents, which would react with MB dye.
- 3) Temperature
- 4) Cation Exchange Capacity (CEC):  $\frac{3 \text{ ml}}{5 \text{ ml}} = 0.6$

# INDIAN INSTITUTE OF TECHNOLOGY

DATE 27/03/25

Exp. 10

Sand Content Kit

SHEET NO. 37

Aim: To measure the sand content in a given trilling fluid using sand content kit.

## Procedure

- 1) Fill the glass measuring tube up to "mud line" with mud
- 2) Use wash bottle to add water to next scribed mark
- 3) Place thumb over mouth of glass measuring tube and shake vigorously
- 4) Pour mixture onto the clean 200 mesh screen.
- 5) Add more water into the glass measuring tube and shake
- 6) Again pour onto 200 mesh screen
- 7) Repeat this procedure until glass measuring tube is clean
- 8) After that wash with water the retained sand present at top of 200 mesh screen.
- 9) Place funnel on top of screen assembly
- 10) Slowly invert assembly and insert tip of funnel into glass measuring tube mouth.
- 11) Spraying a fine spray of fresh water at the back side of the 200 mesh screen that retained sand back into the glass measuring tube mix with water.
- 12) Allow clean sand to settle.
- 13) Read the record vol. Percent sand from graduations on the glass measuring tube.

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

SHEET NO. 38

## Observation and Result

Sand content was found to be = 0.01

## Assignment

- 1) Excessive sand content may harm: All of the above
- 2) Top most line-marking of glass measuring tube: Midline with line ~~no~~
- 3) Water is used for cleaning purpose only if in case of water based mud
- 4) All of the above

Affid,  
02/04/25

(9  
10)

# INDIAN INSTITUTE OF TECHNOLOGY

DATE

Exp10: Retort Kit

SHEET NO.

Aim: To determine the volumes of water, oil and solids in drilling fluid using Retort kit.

Working Principle : The Retort Kit

- This instrument is used to separate water, oil and solid components in drilling fluid samples, and determine their volume percentage of the drilling fluid
- It heats a known vol. of drilling fluid sample in the sample cup, evaporating the liquid components in a distiller, condenser to cool it, collecting in a graduated measuring cylinder.
- Obtain the contents of water, oil and solid components in the drilling fluid sample

Components of Retort Kit :

- Graduated cylinder, 50ml
- Spatula
- Corkscrew
- Square Bar retort Wrench
- Pipe Cleaners

- Steel wool
- Drill Twist
- Brushes - wire and stainless steel
- High Temp. lubricant
- Wetting Agent.

Working Procedure : Retort Kit

- clean and dry the retort assembly and condenser. Condenser passage should be cleaned using a pipe cleaner
- mix the mud sample thoroughly to ensure homogeneity, being careful not to entrain any air or gas entrapment will cause erroneously high retort solids due to reduced mud sample volume.

- Use a clean syringe to take a sample of the mud to be tested
- Fill retort cup slowly to avoid air entrapment. Lightly tap side of cup to expel air. Place lid on the cup. Rotate lid to obtain a proper fit. Be sure a small amount of mud flows out of the hole in the lid. Wipe off any excess mud without wiping any of the samples from inside the cup. Wipe the cup clean with rough cloth.
- Pack retorts body with steel wool.
- Place the retort assembly inside the retort cup. Apply Never-Seal (high temperature lubricant) to the threads of the retort cup with lid in place, hand tighten retort cup onto retort body.
- Apply never-seal to threads of stem and attach condenser.
- Plug cord into 110-volt outlet. Allow the retort to run for a min<sup>m</sup> of 45 mins.
- Unplug retort and allow cooling of retort cup to reach less than 100°F. Disassemble and clean retort assembly and condenser.

Observation

Total liq. volume in the receiver = 35 ml

oil volume = 0 ml

water volume = 35 ml

Calculation

$$1) \text{ Vol. of retort cup, } V_{rc} = \frac{\text{Weight of water (g)}}{\text{Density of water (g/cm}^3\text{)}} = 50 \text{ ml}$$

~~-35 ml.~~

$$2) \text{ Vol. Percent of water, } V_w = \frac{100 \times 35}{50} = 70\%.$$

$$3) \text{ Vol. Percent of oil, } V_o = \frac{100 \times 0}{V_{rc}} = 0\%.$$

$$4) \text{ Vol. Percent of solids, } V_s = \frac{V_{rc}}{100} - (V_w + V_o) = 30\%.$$

Now,  $V_o + V_{Br} + V_{Lars} + V_{wt} = 1$

$$\Rightarrow V_{Br} + V_{Lars} + V_{wt} = 1$$

$$P_{Br} V_o + P_{Br} V_{Br} + P_{Br} V_{Lars} + P_{Br} V_{wt} = MW = 15 \text{ (Given)}$$

Given, MW = 15 PPG

Salinity = 150,000 ppm (15%)

Salt table,  $P_{Br} = 1.116 \text{ g/cc.}$

$$\frac{V_{Br}}{V_w} = \frac{100 / P_{Br}}{100 - \text{wt}\%}$$

$$\therefore V_{Br} = 0.7 \times \left( \frac{100 / 1.116}{100 - 15} \right) = 0.738$$

$$\therefore V_{Lars} + V_{wt} = 1 - 0.738 = 0.262 - (1)$$

$$\text{And, } \rho_{Br} V_{Br} + \rho_{LGS} V_{LGS} + \rho_{wt} V_{wt} = MW = 15$$

$$\Rightarrow 1.116 \times 0.738 + \rho_{LGS} 2.7 V_{LGS} + 4.2 V_{wt} = \frac{15}{8.345} - \textcircled{2}$$

Solving (1) and (2) —

$$V_{LGS} = 0.084 \quad V_{wt} = 0.178$$

~~$$ASGr = \frac{\rho_{LGS} V_{LGS} + \rho_{wt} V_{wt}}{V_{LGS} + V_{wt}} = 3.714$$~~

#### RESULTS:

$$\rightarrow \text{Vol. percent of brine} = 0.738 = 73.8\%$$

$$\rightarrow \text{Vol. percent of oil} = 0\%$$

$$\rightarrow \text{Vol. percent of LGS} = 8.4\%$$

$$\rightarrow \text{Vol. percent of HGS} = 17.8\%$$

$$\rightarrow \text{ASGr of given sample} = 3.714$$

#### PRECAUTIONS:

- The assembly should be cooled to a temp. less than 100°F (37.8°C) from any previous usage
- The threads on the retort should be visually inspected before use of any sign of damage
- The steel wool should be changed out after every test to prevent solids from building up.
- Regular retort should ~~be~~ NEVER be used for drilling fluids containing formate brines

Assignment

Group C

1) Never-seet

3) ~~True~~ False

2) Average specific gravity

4) Retort sample cup

5) Condenser : cools water and oil vapors into liq. for measurement.

Pore cleaner: removes residue from the condenser passage to ensure unobstructed vapor flow.

NumericalWBM  $\rightarrow$  16 lbm/gal having NaCl brineBrine salinity  $\rightarrow$  150000 ppm NaCl (15 wt%)

Retort water content: 33 ml (66% of 50 ml sample)

$$\rho_{\text{brine}} = 2.7 \quad \rho_{\text{barite}} = 4.2$$

$$\rho_{\text{br}} = 1.116 \text{ g/cm}^3$$

(from salt table)

$$V_w = \frac{\text{Vol. of water}}{\text{Vol. rc}} = \frac{33}{50} = 0.66$$

$$\frac{V_{\text{br}}}{V_w} = \frac{100 / \rho_{\text{br}}}{100 - \text{wt}\%} = \frac{100 / 1.116}{100 - 15} = 0.696$$

$$\text{Using, } V_o + V_{\text{bar}} + V_{\text{Lars}} + V_w = 1 \quad [V_o = 0 \because \text{WBM}]$$

$$V_{\text{Lars}} + V_w = 0.304 - (\text{i})$$

# INDIAN INSTITUTE OF TECHNOLOGY

DATE \_\_\_\_\_

SHEET NO. \_\_\_\_\_

$$\rho_0 \cancel{\times}^0 + \rho_{br} V_{br} + \rho_{LGS} V_{LGS} + \rho_{wt} V_{wt} = MW$$

$$\Rightarrow 1.116 \times 0.696 + 2.7 \times V_{LGS} + 4.2 V_{wt} = \frac{16}{8.345} \rightarrow ②$$

Solving ① and ② —

$$V_{LGS} = 0.0908 \quad V_{wt} = 0.2132$$

$$\therefore V_{LGS} = 0.0908 \quad V_{LGS} = V_{wt} = 0.2132 \quad V_{brine} = 0.696$$